



EMINENT  
— SPINE —

# Scoliosis Deformity Pedicle Screw System

Surgical Technique



EMINENT  
— SPINE —

Eminent Spine



EMINENT  
— SPINE —

# Scoliosis Deformity Pedicle Screw System

## Surgical Technique



### **PRODUCT FEATURES**

---

- Aggressive teeth to resist migration
- Sharp Tip to easily penetrate bone
- Central opening for K-wire insertion
- Accepts 1.40mm K-Wire

### **SCREW OPTIONS**

---

#### **PEDICLE SCREW OPTIONS**

| <u>Diameter</u> | <u>Length</u> |
|-----------------|---------------|
| 6.50mm          | 40-60mm   5mm |
| 7.00mm          | 40-60mm   5mm |
| 7.50mm          | 40-60mm   5mm |
| 8.00mm          | 40-60mm   5mm |

#### **ILIAC SCREW OPTIONS**

| <u>Diameter</u> | <u>Length</u>  |
|-----------------|----------------|
| 5.50mm          | 40-110mm   5mm |
| 6.50mm          | 40-110mm   5mm |
| 7.50mm          | 40-110mm   5mm |
| 8.50mm          | 40-110mm   5mm |
| 9.00mm          | 40-110mm   5mm |
| 9.50mm          | 40-110mm   5mm |
| 10.50mm         | 40-110mm   5mm |
| 11.50mm         | 40-110mm   5mm |



# Scoliosis Deformity Pedicle Screw System Surgical Technique



## PRODUCT FEATURES



Closed Side by Side Connector



Open Side by Side Connector



End to End Connector



Lateral Connector

Length: 15-60mm, increments of 5



4.0mm Hex set caps



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 1: Patient Positioning

Place the patient in the prone position on the operating table with hips flexed and the legs adjusted to provide the desired sagittal alignment. Surgical draping and exposure are performed in the routine fashion. Expose the Spine.

Confirm there is enough clearance for a fluoroscopic C-arm to rotate freely for AP, oblique, and lateral views.

## Step 2: Pedicle Identification

The junction of the transverse process and superior articular facet is the typical landmark for entry into the pedicle. Anatomic variations in the individual patients may cause a slight difference of the entry site. Open or MIS options are available to insert the screw.

After the pedicle entry point has been identified, decorticate the entry point with a rongeur or burr, then use a sharp-tipped awl to penetrate the entry point (Figure 1, 2).



Figure 1



Figure 2



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 2: Pedicle Identification

### MODULAR KOCHER

Radiolucent Tip

Small Hole = K-Wire

Middle Hole = Inner Shaft of Targeting Device

Large Hole = Outer Shaft on Targeting Device



Grips the rings on the targeting device,  
does not slip

9 inches long - Gets surgeons hand out of  
X-ray beam

1.4mm K-Wire



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 2: Pedicle Identification

Insert the MIS targeting device through the incision and dock the tip on the pedicle entry point (Figure 3). Using both AP and lateral fluoroscopy, confirm that the appropriate pedicle starting point has been determined. Tap the MIS Targeting Device gently to engage the trocar tip in the pedicle. Insert the trocar tip into the vertebral body. Ensure that the trocar tip does not penetrate the pedicle wall.

Remove the MIS Targeting Device inner stylet (Figure 3, 4).



Figure 3



Figure 4

### MIS TARGETING DEVICE



Inner and Outer Shaft Combined

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 3: Pedicle Probe

Enter the pedicle with a straight or curved probe. Fluoroscopy may be utilized to confirm proper positioning in the sagittal and axial planes (Figures 5, 6).




Figure 5



Figure 6

## LENKE PROBE (STRAIGHT & CURVED)



## Step 4: Ball Tipped Probe

A ball-tipped probe (Depth Sounder) is placed into the pedicle to check for intact walls (Figure 7, 8).



Figure 7



Figure 8

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 5: K-wire Insertion

Insert the K-wire (Figure 9) through the MIS Targeting Device or simply place it in the pedicle hole. Confirm placement with AP and lateral fluoroscopy to ensure that the K-wire does not breach the pedicle or vertebral body wall.



Figure 9



Figure 10

Once the K-wire has been placed to the desired depth, carefully remove the MIS Targeting Device while holding the K-wire in place (Figure 10).

1.4mm Diameter

Blunted or Sharp ends

Stainless Steel

K-wire

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 6: Tap

Insert the Tap over the K-wire (Figure 11). Tap the pedicle with the appropriate diameter tap, which coincides with the screw diameter (Figure 12, 13).



Figure 11




Figure 12



Figure 13

### TAPS



### HANDLES



Straight Handle



T-Handle



Palm Handle

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 7: Screw Selection & Loading

Options include 6.5, 7.0, 7.5 and 8.0mm screw diameter. Lengths range from 40 to 60mm in 5mm increments. The screw is placed on the driver (Figure 14) and the driver is attached to a T-handle or straight ratchet.



### SCREW INSERTERS



Self Retaining



Universal



Universal Removal (3.5 Hex)



Figure 14

## 3.5mm Hex Tip for Screw Insertion

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 8: Screw Insertion

Turn the driver clockwise to insert the screw into the pedicle. The screw should be parallel to the endplates and extend 50-80% into the vertebral body. Fluoroscopy should be utilized to confirm proper positioning (Figure 15, 16).




1st Screw



Figure 15

Repeat same process for the additional screws.



2nd Screw



Figure 16

## 3.5mm Hex Tip for Screw Insertion



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 9: Rod Selection & Placement

5.5mm straight and pre-bent rods are available in lengths starting at 40mm.

### Rod

Bullet nose for  
easier passage



Visible line marking for lordosis

### ROD HOLDERS

6 Options:

Rod Holders (Open)



Rod Holders (MIS)



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 10: Rod Selection & Placement

After all the pedicle screws are in position, insert the appropriate rod clamped by a rod holder and seat into the screw tulip (Figure 17).




Figure 17

A rod bender is available in needed (Figure 18).



Figure 18

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 11: Pedicle Set Screw Caps



4.0mm Hex Tip for Set Screw

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 11 Cont: Pedicle Set Screw Caps

Six options (all with a tapered tip) are available to insert the Pedicle Set Screw Caps.



### 6 Options:



Palm Handled Welded Together



Straight Handled Welded Together



### 4.0mm Hex Tip for Set Screw

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 11 Cont: Pedicle Set Screw Caps

Set cap screws are inserted into the threaded portion of the screw tulip using the set cap starter (Figure 19). The set cap screws should not be fully tightened but should remain loose, so that the rod can be placed in the appropriate position. After the rod is in the correct position, one of the set screws should be provisionally tightened (Figure 20).



Figure 19



Figure 20

## Step 12: Compression/Distrraction

If compression/ distraction is needed, provisionally tighten one set screw, and slightly loosen the other set screw to allow for the rod to slide. The compression/distraction is accomplished to the correct frontal and/or sagittal plane, then the set screw is provisionally tightened (Figure 21, 22).




Figure 21



Figure 22



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 13: Final Tightening

Counter Torque and Torque Wrench are available for final tightening of the pedicle set screw caps.

### COUNTER TORQUE



MIS & OPEN  
Visible markings for scrub  
tech to easily identify



### TORQUE WRENCH



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 13 Cont: Final Tightening

The locking set screws are tightened to the rod using the counter-torque wrench, which is seated on the rod. The driver is then inserted through the counter-torque wrench onto the set screws. The final tightening device is a torque-limiting wrench preset to 10Nm. Turn the torque handle clockwise until a click is heard. The set screw is then fully tightened. Repeat on the remaining screws (Figure 23).



Figure 23



# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 14: Final Position

Obtain Anterior, Oblique and Lateral x-rays after final position is complete (Figure 24).



Figure 24

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Step 15: Transverse Connector

Transverse connectors may be placed to connect the rods together. They may be placed between the rods to enhance torsional stability of the overall construct. Three different lengths are available. The middle screw should be loosened to allow for length adjustment. The outer two screws should be loosened to the point that the end of the connector snaps over the rod. Tighten the outer 2 screws with a torque limiting wrench (preset to 3.5N/m), and then tighten the middle screw. Turn the wrench clockwise until a click is heard (Figure 25, 26).



Figure 25




Figure 26

## Step 16: Closure

Wound closure is then performed in the customary manner.

## Step 17: Revision/Removal Procedure

Remove the transverse connector, followed by the set screws by turning the set screw driver counter clockwise, then remove the rod. Insert the pedicle screw driver into the pedicle screw and turn the driver counter clockwise and remove the pedicle screw. If revision is required, a larger diameter pedicle screw may be needed.

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Closed Side by Side Connector



Attach Closed Side by Side Connector to Rod. Tighten Set Screw Cap on the side connected to the Rod (Figure 1).



Figure 1



Insert Rod to Closed Side by Side Connector and Pedicle Screw (Figure 2) and tighten remaining Set Screw Cap to secure (Figure 3).



Figure 2

Figure 3

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Closed Side by Side Connector



Figure 4

Insert 4.0mm Hex Set Screw Cap  
on opposite end of the Rod to  
secure (Figure 4).

Tighten the 4.0mm Hex to  
10N-m (Figure 5).



Figure 5

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Open Side by Side Connector



Figure 1



Attach Open Side by Side Connector to Rod (Figure 1).



Figure 2

Insert Rod and tighten the 4.0mm Hex to 5N-m (Figure 2).

# Scoliosis Deformity Pedicle Screw System Surgical Technique

## End to End Connector



Figure 1



Attach Connector to Rod (Figure 1).



Insert Rod to Connector and Pedicle Screw (Figure 2).



Figure 2



Figure 3

Insert Rod and tighten the 4.0mm Hex to 3.5N-m (Figure 3).

# Scoliosis Deformity Pedicle Screw System Surgical Technique



Lateral Connector

Iliac Fixation:

Iliac screw entry point is medial to the medial iliac crest wall, 2cm below PSIS, and lateral to S1 foramen (Figure 1). The trajectory is 45° caudally and laterally towards the anterior inferior iliac spine (Figure 2).



Figure 1



Figure 2

It is recommended to notch the iliac crest sufficiently to recess the tulip, to avoid implant prominence (Figure 3).



Figure 3

# Scoliosis Deformity Pedicle Screw System Surgical Technique



Insert the Lateral Connector to the Iliac Screw and secure with 4.0mm Hex Set Screw Cap (Figure 4).



Figure 4

Insert Rod into Pedicle Screw Tulip and Lateral Connector Tulip and provisionally tighten with 4.0mm Hex Set Screw Caps (Figure 5).



Figure 5

Once the set caps are tightened fully, remove the set cap driver (Figure 6).



Figure 6



# Scoliosis Deformity Pedicle Screw System Surgical Technique



Lateral Connector



Tighten the 4.0mm Hex to 10N-m utilizing the Counter Torque on the pedicle screws (Figure 7).

Figure 7

Tighten the Iliac Rod Connector to 10N-m utilizing the Counter Torque (Figure 8).



Figure 8

# Scoliosis Deformity Pedicle Screw System Surgical Technique



## Step : Closure

Wound closure is then performed in the customary manner.



## Step : Revision/Removal Procedure

Remove the set screw caps by turning the driver counter clockwise, then remove the rod connectors, rods and/or connectors. Insert the pedicle screw driver into the screw and turn the driver counter clockwise and remove the screw. If revision is required, a larger diameter screw may be needed.





# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Pedicle System Instrument Tray





# Scoliosis Deformity Pedicle Screw System Surgical Technique

## Pedicle System Screw Tray



# Scoliosis Deformity Pedicle Screw System Instrument Catalog

| Part Number | Pedicle 3D Screw System      | Materials  | Notes        | Image   |
|-------------|------------------------------|--|--------------|---|
| PS-105      | Pedicle Instrument Tray      | Anodized Aluminum base, includes Nylon coated SS brackets, anodized aluminum post, electropolished TRI-corners, and 1-color silk-screened graphics on bottom and side panels | 3 Level Tray |    |
| PI-106      | Pedicle Open Implant Tray    | Anodized Aluminum base, includes Nylon coated SS brackets, anodized aluminum post, electropolished TRI-corners, and 1-color silk-screened graphics on bottom and side panels | 2 Level Tray |    |
| DB-1003L    | Pedicle Screw Driver, Long   | 316 Stainless Steel per ASTM A276  |              |    |
| DB-1006     | Set Cap Caddy                | T6061 Aluminum, Smoke Polyphenylsulfone Sheet (Radel R5000)  |              |    |
| DB-1009     | AO Cross Link Driver         | 455 Stainless Alloy per ASTM F899  |              |  |
| DB-1010     | Cross Link Caddy             | T6061 Aluminum, Smoke Polyphenylsulfone Sheet (Radel R5000)  |              |  |
| DB-1012     | Set Cap Starter, Double Flex | 455 Stainless Alloy per ASTM F899  |              |  |
| DB-1014     | Radiolucent Modular Kocher   | Medical Grade Polypropylene (Tecapo MT)  |              |  |
| DB-1020     | Straight Bone Probe          | 17-4 Stainless Steel per ASTM A564   |              |  |
| DB-1021     | Curved Bone Probe            | 17-4 Stainless Steel per ASTM A564   |              |  |
| DB-1022     | Pedicle Awl                  | 17-4 Stainless Steel per ASTM A564   |              |  |
| DB-1023     | Counter Torque - Open        | 17-4 Stainless Steel per ASTM A564   |              |  |
| DB-1023L    | Counter Torque - Open Long   | 17-4 Stainless Steel per ASTM A564   |              |  |

|          |  |                                    |   |
|----------|--|------------------------------------|---|
| DB-1025  | Wall Check, Depth Sounder                | 17-4 Stainless Steel per ASTM A564 |    |
| DB-1031  | Counter Torque, MIS                      | 17-4 Stainless Steel per ASTM A564 |    |
| DB-1031L | Counter Torque, MIS Long                 | 17-4 Stainless Steel per ASTM A564 |    |
| DB-1033  | Rod Inserter, MIS                        | 17-4 Stainless Steel per ASTM A564 |    |
| DB-1034L | Pedicle Screw Driver, Threaded Universal | 455 Stainless Alloy per ASTM F899  |    |
| DB-1036  | Pedicle Distractor                       | 17-4 Stainless Steel per ASTM A564 |    |
| DB-1037  | Pedicle Compressor                       | 17-4 Stainless Steel per ASTM A564 |   |
| DB-1050L | Final Cap Driver, 4.0mm Long             | 455 Stainless Alloy per ASTM F899  |  |
| DB-1080L | Pedicle MIS Target Tap System, Long      | 17-4 Stainless Steel per ASTM A564 |  |
| DB-1083L | Set Cap Starter, 4.0mm Flex Tip, Long    | 455 Stainless Alloy per ASTM F899  |  |
| DB-1085  | Rod Bender                               | 17-4 Stainless Steel per ASTM A564 |  |
| DB-1086  | Rod Holder A, Open                       | 17-4 Stainless Steel per ASTM A564 |  |
| DB-1105L | Cap Starter, 4.0mm, Long                 | 455 Stainless Alloy per ASTM F899  |  |
| DB-1108  | Rod Inserter, Version-K                  | 17-4 Stainless Steel per ASTM A564 |  |
| DB-1109  | Pedicle Lenke Probe, Straight            | 17-4 Stainless Steel per ASTM A564 |  |
| DB-1110  | Pedicle Lenke Probe, Curved              | 17-4 Stainless Steel per ASTM A564 |  |

# Scoliosis Deformity Pedicle Screw System Instrument Catalog

|          |  |                                    |  |   |
|----------|--|------------------------------------|--|---|
| DB-1114  | Pedicle Screw Driver - Retaining       | 455 Stainless Alloy per ASTM F899  |  |                                        |
| DB-1145  | 4.5mm Tap, Cannulated                  | 17-4 Stainless Steel per ASTM A564 |  |                                        |
| DB-1155  | 5.5mm Tap, Cannulated                  | 17-4 Stainless Steel per ASTM A564 |  |   |
| DB-1165  | 6.5mm Tap, Cannulated                  | 17-4 Stainless Steel per ASTM A564 |  |   |
| DB-1175  | 7.5mm Tap, Cannulated                  | 17-4 Stainless Steel per ASTM A564 |  |   |
| ES-802   | 1/4 Drive Ratchet, T-Handle            | 17-4 Stainless Steel per ASTM A564 |  |   |
| ES-803   | 1/4 Drive Ratchet, I-Handle            | 17-4 Stainless Steel per ASTM A564 |  |                                        |
| ES-804   | 1/4 Drive Ratchet, Palm Handle         | 17-4 Stainless Steel per ASTM A564 |  |                                        |
| ES-805   | 1/4 Drive 10Nm Final Torque            | 17-4 Stainless Steel per ASTM A564 |  |                                      |
| ES-806   | A/O Drive 3.5NM Torque                 | 17-4 Stainless Steel per ASTM A564 |  |                                      |
| ESP-0001 | K-Wire, 1.4mm x 18", 316 SS, Flat Ends | 316 Stainless Steel per ASTM A276  |  |  |